

**THE PILL TIMER™ IS A MEDICATION MANAGEMENT SYSTEM DESIGNED TO HELP OUR PATIENTS  
RECEIVE THE FULL THERAPEUTIC BENEFIT FROM THEIR PRESCRIPTION MEDICATION  
WHILE REDUCING THE RISK OF TREATMENT FAILURE AND ADVERSE DRUG EVENTS**

**Simple to Set - Easy to Use**

1. Remove Pill Timer from Prescription Vial.



2. Press and Hold button until prescribed dosage interval is displayed on screen.



3. Place Pill Timer securely on Prescription Vial.



**Timer automatically resets to 0:00 after each use.**

**To Test The Pill Timer™**

1. Place Pill Timer securely on prescription vial.
2. Press and HOLD Green button for 8 Seconds.

**Prescribed Medication Regimen**

<b>Once Daily:</b>	Set to <b>24</b>
<b>Twice Daily:</b>	Set to <b>12</b>
<b>Three Times per day:</b>	Set to <b>8</b>
<b>Four Times per day:</b>	Set to <b>6</b>
<b>As Needed (No Alert):</b>	Set to <b>0</b>

- The Pill Timer is operating if 2 dots ( : ) between HRS : MIN are blinking when properly placed on prescription vial.
- If cap is not removed from vial for 7 days; LCD will go blank and Pill Timer will go into sleep mode. To reactivate Pill Timer simply press button or remove from vial.
- Battery is designed to operate for life of prescription. Please discard Pill Timer after 90 days or if battery is no longer functioning.
- The Pill Timer™ should be discarded if LCD is blank, displays 88:88, or 24.

- **When Pill Timer is on vial**, press button to display hourly dosage interval. This will not change setting or stop Pill Timer.
- **Counter Feature:** Hold button for 4 seconds to display how many times Pill Timer has been removed from vial.
- **Reminder Feature:** The Pill Timer™ has a built in hourly reminder alarm that will activate every hour dosage has been missed; up to 24 hrs.
- **Snooze Feature:** Press Green Button to stop reminder alert. This will not affect subsequent reminder alarms.

The Pill Timer™ has a 90 day manufacturer's warranty. This warranty covers defects or malfunctions only. The warranty does not apply to battery life, loss or damage to the product due to abuse, misuse or exposure to extreme temperatures, water or any other liquids. There are no user serviceable parts. Do not attempt to take The Pill Timer™ apart. The warranty will be voided for any of these instances. The manufacturer's warranty is limited solely to the replacement of any defective unit. Defective units should be returned with a dated sales receipt for replacement to: Med Time Technology, Inc. 6464 NW 5th Way Fort Lauderdale, FL 33309. (561.) 988-1888. If The Pill Timer™ is defective, Please Return To Your Pharmacy For A Free Replacement.